
ASPEN[®]CTO

CERVICAL THORACIC ORTHOSIS
ZERVIKAL-THORAKALE ORTHESE
ÓRTESIS TORÁCICA CERVICAL
ORTESI CERVICO-TORACICA
ORTHÈSE CERVICO-THORACIQUE
CERVICALE THORACALE ORTHÈSE
KRČNI HRUDNÍ ORTÉZA
ЦЕРВИКАЛ ТHOPAKИH OPTHЭZ
CERVİKAL BRYSTÖRTROSE
ORTEZA SZYJNO-PIERSIOWA
KAULARANKA-RINKEHÄORTOSI
CERVİKAL THORAKAL ORTOS
ОРТЕЗ ДЛЯ ШЕЙНО-ГРУДНОГО ОТДЕЛА ПОЗВОНОЧНИКА
ÓRTESE TORÁCICA CERVICAL
頸椎矫形器

ENGLISH
1. SIZING
 A. Use the Aspen sizing guide to measure the patient for the correct size Aspen Collar. The sizing guide can be found at aspennp.com.
 B. Meet de borstomtrek 2,5 cm onder het uiteinde van de Xiphoid process.
 C. For circumferences smaller than 39", trim the ends of the belt panels at the notch.
 D. Attach the belt by aligning the keyhole over the aluminum rivet and pull to lock in place. Do this on both sides.
2. FITTING
 A. Fit the collar to the patient as described in the Aspen collar instruction sheet. The instruction sheet can be found at aspennp.com.
 B. (Supine): Side the CTO under the patient while compressing the mattress cover. Layroll the patient if needed.) Line up the control clamp to the centerline on the chest and attach the hook and loop.
 C. Snap the shoulder straps into place and loosely tighten. Bars may be bent to accommodate patient anatomy.
 NOTE: To relieve pressure at the front on the collar, bend bars away from the patient.
3. CONTROL CLAMP
 Use a phillips head screw driver to loosen screws to adjust the height of the bars. Tighten the clamp after adjustment.
4. ADJUSTMENT
 A. Slide the support posts on the front of the CTO upward until the patient's head is properly supported. Tighten the screws on the support clamp.
 B. Adjust the collar back in the same manner. (Hint: the higher the back of the collar, the more effective the motion restriction.)
 C. Check for symmetry and tighten shoulder straps removing slack.
NOTES
 • The Aspen[®] CTO can be altered to mirror changes in the patient's condition. For supine patients where posterior support posts are desirable, a pillow may be used under one shoulder to minimize pressure caused by the posterior support post.
 • Be aware that excessive tightening of the shoulder straps will cause traction.
 • Do not use a full size screwdriver to tighten or adjust the support post clamps while the device is on the patient.
 • When advancing to a 4-post configuration, ensure that the post strap covers the support post.
 • For kyphotic patients or for others with postural abnormalities, it is possible to bend the support posts for a better fit. To accomplish this, use a pair of bending irons or remove the posts from the clamp and bend using a vise to ensure that they remain symmetrical. Be careful not to weaken the posts by bending them repeatedly or by bending them more than 45 degrees.
 • Replacement pads are available from your doctor, local orthotic and prosthetic facility, or by calling Aspen Medical Products at 800-295-2776.
INTENDED USE
 To provide motion restriction of the cervical spine.
INDICATIONS
 Decompressive procedures, Fracture management, Post-surgical stabilization, Whiplash, Spains/strains, Cervicogenic headache, Cervical disc syndrome, Musculoskeletal pain, Radiculopathy.
CONTRAINDICATIONS
 Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.
PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.
CARE INSTRUCTIONS
 Hand wash only at 30° C. Wash with mild detergent; Air dry. Do not use bleach or other chemicals.
MATERIALS
 HDPE, Nylon, Polyurethane, Cotton, Aluminum
NOT MADE WITH NATURAL RUBBER LATEX.
DISPOSAL
 Product may be safely disposed in accordance with local laws.
REPORTING
 Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).
LIMITED WARRANTY
 Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).
 Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
 (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
 (b) The product is altered;
 (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
 (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.
 THE REMEDY UNDER THIS WARRANTY IS LIMITED TO REPAIR OR REPLACEMENT OF THE DEFECTIVE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.
 This warranty gives you specific legal rights and you may have other rights which vary from location to location.
 A listing of authorized distributors in your area is available upon request.

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Dans les cas suivants, Aspen Medical Products, LLC n'aura aucune obligation au titre de cette garantie limitée:
(a) Le produit n'a pas été acheté auprès d'Aspen Medical Products, LLC ou via ses canaux de distribution autorisés ;
(c) Toute pièce non fournie par Aspen Medical Products, LLC a été insérée dans le produit ; ou
(d) Le produit n'est pas utilisé conformément au mode d'emploi d'Aspen Medical Products, LLC.